

 Liquid Media	BML01	BML02	BML03
Amino Acids	mg/L	mg/L	mg/L
L-Arginine, Hydrochloride		2106.0000	21.0600
L-Cystine, Dihydrochloride		1565.0000	15.6500
L-Histidine, Hydrochloride, Monohydrate		1050.0000	10.5000
L-Isoleucine		2623.0000	26.2300
L-Leucine		2623.0000	26.2300
L-Lysine, Hydrochloride		3647.0000	36.4700
L-Methionine		746.0000	7.4600
L-Phenylalanine		1651.0000	16.5100
L-Threonine		2382.0000	23.8200
L-Tryptophan		400.0000	4.0000
L-Tyrosine, Disodium, Dihydrate		2595.0000	25.9500
L-Valine		2343.0000	23.4300
Inorganic Salts	mg/L	mg/L	mg/L
Calcium Chloride, Anhydrous			200.0000
Magnesium Sulfate, Anhydrous			97.6700
Potassium Chloride			400.0000
Sodium Chloride			6800.0000
Sodium Phosphate, Monobasic, Monohydrate			140.0000
Other Components	mg/L	mg/L	mg/L
D-Glucose, Anhydrous			1000.0000
HEPES			4770.0000
Phenol Red, Sodium Salt			17.0000



Liquid Media

BML01

BML02

BML03

Vitamins

mg/L

mg/L

mg/L

Biotin	100.0000		1.0000
Choline Chloride	100.0000		1.0000
D-Calcium Pantothenate	100.0000		1.0000
Folic Acid	100.0000		1.0000
Myo-Inositol	200.0000		2.0000
Niacinamide	100.0000		1.0000
Pyridoxal, Hydrochloride	100.0000		1.0000
Riboflavin	10.0000		0.1000
Thiamine, Hydrochloride	100.0000		1.0000