




# Comparative Formulation

DKW

 Powdered Media	DKP01	DKP02
<b>Amino Acids</b>	mg/L	mg/L
Glycine		2.0000
<b>Macronutrients</b>	mg/L	mg/L
Ammonium Nitrate	1416.0000	1416.0000
Calcium Chloride, Anhydrous	112.5000	112.5000
Calcium Nitrate, Tetrahydrate	1367.0000	1367.0000
Magnesium Sulfate, Anhydrous	361.4900	361.4900
Potassium Phosphate, Monobasic, Anhydrous	265.0000	265.0000
Potassium Sulfate, Anhydrous	1559.0000	1559.0000
<b>Micronutrients</b>	mg/L	mg/L
Boric Acid	4.8000	4.8000
Cupric Sulfate, Pentahydrate	0.2500	0.2500
EDTA, Disodium Salt, Dihydrate	45.4000	45.4000
Ferrous Sulfate, Heptahydrate	33.8000	33.8000
Manganese Sulfate, Monohydrate	33.5000	33.5000
Molybdic Acid Sodium Salt, Dihydrate	0.3900	0.3900
Nickel Sulfate, Hexahydrate	0.0050	0.0050
Zinc Nitrate, Hexahydrate	17.0000	17.0000



## Powdered Media

DKP01

DKP02

### Vitamins

mg/L

mg/L

Myo-Inositol

100.0000

Nicotinic Acid

1.0000

Thiamine, Hydrochloride

2.0000