

## Earle's Balanced Salts

### Product Number EBP01-10LT

With calcium, magnesium, and phenol red. Without sodium bicarbonate.

| Components                                                                                      | mg/L      |
|-------------------------------------------------------------------------------------------------|-----------|
| Calcium Chloride, Dihydrate ( $\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$ )                       | 264.9000  |
| D-Glucose, Anhydrous ( $\text{C}_6\text{H}_{12}\text{O}_6$ )                                    | 1000.0000 |
| Magnesium Sulfate, Anhydrous ( $\text{MgSO}_4$ )                                                | 97.7000   |
| Phenol Red, Sodium Salt ( $\text{C}_{19}\text{H}_{13}\text{NaO}_5\text{S}$ )                    | 17.0000   |
| Potassium Chloride (KCl)                                                                        | 400.0000  |
| Sodium Chloride (NaCl)                                                                          | 6800.0000 |
| Sodium Phosphate, Monobasic, Monohydrate ( $\text{NaH}_2\text{PO}_4 \cdot \text{H}_2\text{O}$ ) | 140.0000  |