

Hank's Balanced Salts

Product Number HBL02-6X500ML

HBS solution with calcium, magnesium, and phenol red.

Sterile filtered. Store at room temperature.

Components	mg/L
Calcium Chloride, Anhydrous (CaCl_2)	140.0000
D-Glucose, Anhydrous ($\text{C}_6\text{H}_{12}\text{O}_6$)	1000.0000
Magnesium Sulfate, Anhydrous (MgSO_4)	98.0000
Phenol Red, Sodium Salt ($\text{C}_{19}\text{H}_{13}\text{NaO}_5\text{S}$)	10.0000
Potassium Chloride (KCl)	400.0000
Potassium Phosphate, Monobasic, Anhydrous (KH_2PO_4)	60.0000
Sodium Chloride (NaCl)	8000.0000
Sodium Phosphate Dibasic, Anhydrous (Na_2HPO_4)	48.0000