


 Powdered Media	HSP01	HSP02	HSP03
Amino Acids	mg/L	mg/L	mg/L
Glycine	2.0000	2.0000	
Macronutrients	mg/L	mg/L	mg/L
Ammonium Nitrate	1650.0000	1650.0000	1650.0000
Calcium Chloride, Anhydrous	332.2000	332.2000	332.2000
Magnesium Sulfate, Anhydrous	180.7000	180.7000	180.7000
Potassium Nitrate	1900.0000	1900.0000	1900.0000
Potassium Phosphate, Monobasic, Anhydrous	300.0000	300.0000	300.0000
Sodium Phosphate, Monobasic, Monohydrate	170.0000	170.0000	170.0000
Micronutrients	mg/L	mg/L	mg/L
Boric Acid	6.2000	6.2000	6.2000
Cobalt Chloride, Hexahydrate	0.0250	0.0250	0.0250
Cupric Sulfate, Pentahydrate	0.0250	0.0250	0.0250
EDTA, Disodium Salt, Dihydrate	37.2600	37.2600	37.2600
Ferrous Sulfate, Heptahydrate	27.8000	27.8000	27.8000
Manganese Sulfate, Monohydrate	16.9000	16.9000	16.9000
Molybdic Acid Sodium Salt, Dihydrate	0.2500	0.2500	0.2500
Potassium Iodide	0.8300	0.8300	0.8300
Zinc Sulfate, Heptahydrate	8.6000	8.6000	8.6000

 Powdered Media	HSP01	HSP02	HSP03
Other Components	mg/L	mg/L	mg/L
6-Benzylaminopurine	2.0000	0.1000	0.1000
Adenine Sulfate, Dihydrate	160.0000	160.0000	
Casein Digest Peptone (Tryptone)	500.0000	500.0000	500.0000
Naphthaleneacetic Acid	0.5000	0.5000	0.5000
Vitamins	mg/L	mg/L	mg/L
Myo-Inositol	100.0000	100.0000	100.0000
Thiamine, Hydrochloride	0.4000	0.4000	0.4000