

 Liquid Media	ITL01	ITL02	ITL03	ITL04
Inorganic Salts	mg/L	mg/L	mg/L	mg/L
Calcium Chloride, Anhydrous		200.0000	200.0000	200.0000
Magnesium Sulfate, Heptahydrate		200.0000	200.0000	200.0000
Potassium Chloride	400.0000	400.0000	400.0000	400.0000
Sodium Bicarbonate	2200.0000	2200.0000	2200.0000	2200.0000
Sodium Chloride	6800.0000	6800.0000	6800.0000	6800.0000
Sodium Phosphate, Monobasic, Monohydrate	140.0000	140.0000	140.0000	140.0000
Sodium Selenite, Anhydrous	0.6700	0.6700	0.6700	0.6700
Other Components	mg/L	mg/L	mg/L	mg/L
D-Glucose, Anhydrous	1000.0000	1000.0000	1000.0000	1000.0000
Ethanolamine		200.0000		
Insulin, Recombinant Human	1000.0000	1000.0000	1000.0000	1000.0000
Pyruvic Acid, Sodium Salt	11000.0000			11000.0000
Transferrin	550.0000	550.0000	550.0000	550.0000