



CERTIFICATE OF ANALYSIS

Product Name: KSOM+AA
Catalog Number: IVL04 Lot Number: 01161007

Test	Specifications	Results
Appearance	Clear, Colorless Solution	PASS
Sterility		
Aerobic	Negative	NEGATIVE
Anaerobic	Negative	NEGATIVE
pH <i>(pH was measured at ~ 22° C).</i>	7.30 - 7.40	7.40
Osmolality	Record mM/Kg	262 mM/Kg
Endotoxin	≤ 10.0 EU/mL <i>Endotoxin testing was performed using the LAL detection method.</i>	< 0.1 EU/mL
Mycoplasma	Negative <i>Mycoplasma test was performed using MycoAlert® Mycoplasma Detection Assay.</i>	NEGATIVE

Manufacture Date: Jan 2016
Expiration Date: 31 Jan 2018

This product was produced and tested in accordance with all relevant SOP's and meets or exceeds the quality standards established by Caisson Laboratories, Inc.

This product is intended for laboratory use only.

Kurt Rieben 1-28-16
Quality Control