



www.caissonlabs.com

CERTIFICATE OF ANALYSIS

KSOM+AA

KSOM with Amino Acids

Catalog Number: IVL04

Lot Number: 01211011

MADE IN THE USA

Storage: 2° to 8°C

Manufacture Date: Jan-2021

Expiration Date: **31-Jul-2022**

Test	Specifications	Results
Appearance	Clear, Colorless Solution	Pass
Sterility (<i>USP <71></i>)		
Aerobes, Anaerobes,	Negative	Negative
Fungi		
pH (<i>measured at 25°C</i>)	7.3 – 7.4	7.4
Osmolality	Record mM/Kg	259 mM/Kg
Endotoxin	≤ 10.0 EU/mL	< 0.1 EU/ml
<i>Endotoxin Testing was Performed using the LAL Detection Method</i>		
Mycoplasma	Negative	Negative
<i>Mycoplasma Test was Performed using MycoAlert® Mycoplasma Detection Assay</i>		

This product was produced and tested in accordance with all relevant SOP's and meets or exceeds the quality standards established by Caisson Laboratories, Inc.

This product is intended for laboratory use only.

Prepared by:

Kurt Rieben
QC Technician

Reviewed by:

BSE/TSE Declaration

To the present state of our knowledge we declare that this product does not contain any BSE/TSE AGENT. This declaration is given in good faith and no warranty express or implied with respect to quality and properties is made.