



CERTIFICATE OF ANALYSIS

Product Name: KSOM+AA

Catalog Number: IVL04

Lot Number: 11171007

MADE IN THE USA

Storage: 2 to 8°C

Manufacture Date: NOV-2017

Expiration Date: 30-NOV-2019

Test	Specifications	Results
Appearance	Clear, Colorless Solution	PASS
Sterility		
Aerobic	Negative	NEGATIVE
Anaerobic	Negative	NEGATIVE
pH (measured at 22°C)	7.30 - 7.40	7.39
Osmolality	Record mM/Kg	260 mM/Kg
Endotoxin	≤ 10.0 EU/mL	< 1 EU/mL
<i>Endotoxin testing was performed using the LAL detection method.</i>		
Mycoplasma	Negative	NEGATIVE
<i>Mycoplasma test was performed using MycoAlert® Mycoplasma Detection Assay.</i>		

This product was produced and tested in accordance with all relevant SOP's and meets or exceeds the quality standards established by Caisson Laboratories, Inc.

This product is intended for laboratory use only.

Prepared by:

Reviewed by:

Kurt Rieben
Quality Control

27-NOV-2017

BSE/TSE Declaration

To the present state of our knowledge we declare that this product does not contain any BSE/TSE AGENT. This declaration is given in good faith and no warranty express or implied with respect to quality and properties is made.