



CERTIFICATE OF ANALYSIS

Product Name: LB (Lennox) Agar

Catalog Number: LBP04, YQP01

Lot Number: 09200015

Storage: 15 to 30°C

Manufacture Date: SEP-2020

Expiration Date: 30-SEP-2022

MADE IN THE USA

Test	Specifications	Results
Appearance	Off-White to Yellow / Tan Powder	PASS
Solubility (in H ₂ O)	35.00 g/L in ≤ 30 minutes (suspension)	20 minutes
Moisture Content	Record %	7.4%
Particle Size (100 mesh)	≥ 90%	≥ 90%
pH (measured at 28°C)	6.8 - 7.2	6.87
Osmolality	Record mM/Kg	240 mM/Kg
Growth Promotion	Supports growth of bacterial strains: <i>Escherichia coli</i> and <i>Staphylococcus epidermidis</i>	PASS

This product was produced and tested in accordance with all relevant SOP's and meets or exceeds the quality standards established by Caisson Laboratories, Inc.

This product is intended for laboratory use only.

Prepared by:

Kurt Rieben
Quality Control

07-OCT-2020

Reviewed by:

BSE/TSE Declaration

To the present state of our knowledge we declare that this product does not contain any BSE/TSE AGENT. This declaration is given in good faith and no warranty express or implied with respect to quality and properties is made.