



www.caissonlabs.com

CERTIFICATE OF ANALYSIS

MCDB 131 Medium

Without L-Glutamine

Catalog Number: MBL02
Lot Number: 11211001

MADE IN THE USA

Storage: 2° to 8°C
Manufacture Date: Nov-2021
Expiration Date: **30-Nov-2023**

Test	Specifications	Results
Appearance	Clear, Red Solution	Pass
Sterility (<i>USP <71></i>)		
Aerobes, Anaerobes,	Negative	Negative
Fungi		
pH (<i>measured at 25°C</i>)	7.2 – 7.6	7.59
Osmolality	260 – 300 mM/Kg	269 mM/Kg
Endotoxin	≤ 1.0 EU/mL	< 0.1 EU/ml
<i>Endotoxin Testing was Performed using the LAL Detection Method</i>		
Mycoplasma	Negative	Negative
<i>Mycoplasma Test was Performed using MycoAlert® Mycoplasma Detection Assay</i>		
Animal Cell Culture		
Cell Growth Promotion	3 Passages	Pass
Viability	≥ 80 %	≥ 85 %
<i>Medium was supplemented with 10% FBS and 292 mg/L L-Glutamine, then tested with VERO cells grown in a monolayer in a 5% CO₂ atmosphere.</i>		

This product was produced and tested in accordance with all relevant SOP's and meets or exceeds the quality standards established by Caisson Laboratories, Inc.

This product is intended for laboratory use only.

Prepared by:

Alisha Wright
QC Technician

Reviewed by:

15-Nov-2021

BSE/TSE Declaration

To the present state of our knowledge we declare that this product does not contain any BSE/TSE AGENT. This declaration is given in good faith and no warranty express or implied with respect to quality and properties is made.