



CERTIFICATE OF ANALYSIS

Product Name: RPMI 1640 Medium
With L-Glutamine and 25mM HEPES

Catalog Number: RPL09 Lot Number: 03161013

Test	Specifications	Results
Appearance	Clear, Red-Orange Solution	PASS
Sterility		
Aerobic	Negative	NEGATIVE
Anaerobic	Negative	NEGATIVE
pH <i>(pH was measured at ~ 22 ° C).</i>	7.00 - 7.60	7.39
Osmolality	250 - 320 mM/Kg	308 mM/Kg
Endotoxin	≤ 1.0 EU/mL	< 0.1 EU/mL
<i>Endotoxin testing was performed using the LAL detection method.</i>		
Mycoplasma	Negative	NEGATIVE
<i>Mycoplasma test was performed using MycoAlert® Mycoplasma Detection Assay.</i>		
Animal Cell Culture		
Cell Growth Promotion	3 Passages	PASS
Viability	≥ 80%	≥ 85%
<i>Medium was supplemented with 10% FBS, then tested using BHK cells grown in a monolayer in a 5% CO₂ atmosphere.</i>		

Manufacture Date: Mar 2016

Expiration Date: 30 Sep 2017

This product was produced and tested in accordance with all relevant SOP's and meets or exceeds the quality standards established by Caisson Laboratories, Inc.

This product is intended for laboratory use only.

Kurt Rieben
Quality Control

3-24-16