



CERTIFICATE OF ANALYSIS

Product Name: RPMI 1640 Medium
With L-Glutamine and 25mM HEPES

MADE IN THE USA

Catalog Number: RPL09

Storage: 2 to 8°C

Manufacture Date: MAR-2017

Lot Number: 03171004

Expiration Date: 30-SEP-2018

Test	Specifications	Results
Appearance	Clear, Red-Orange Solution	PASS
Sterility		
Aerobic	Negative	NEGATIVE
Anaerobic	Negative	NEGATIVE
pH (measured at ~ 22°C)	7.00 - 7.60	7.24
Osmolality	250 - 320 mM/Kg	280 mM/Kg
Endotoxin	≤ 1.0 EU/mL	< 0.1 EU/mL
<i>Endotoxin testing was performed using the LAL detection method.</i>		
Mycoplasma	Negative	NEGATIVE
<i>Mycoplasma test was performed using MycoAlert® Mycoplasma Detection Assay.</i>		
Animal Cell Culture		
Cell Growth Promotion	3 Passages	PASS
Viability	≥ 80%	> 90%
<i>Medium was supplemented with 10% FBS, then tested using VERO cells grown in a monolayer in a 5% CO₂ atmosphere.</i>		

This product was produced and tested in accordance with all relevant SOP's and meets or exceeds the quality standards established by Caisson Laboratories, Inc.

This product is intended for laboratory use only.

Prepared by:

Reviewed by:

Kurt Rieben
Quality Control

3-23-17