




Comparative Formulation

Chu's

 Powdered Media	CHP01	CHP03
Amino Acids	mg/L	mg/L
Glycine		2.0000
Macronutrients	mg/L	mg/L
Ammonium Sulfate	463.0000	463.0000
Calcium Chloride, Anhydrous	125.3300	125.3300
Magnesium Sulfate, Anhydrous	90.3700	90.3700
Potassium Nitrate	2830.0000	2830.0000
Potassium Phosphate, Monobasic, Anhydrous	400.0000	400.0000
Micronutrients	mg/L	mg/L
Boric Acid	1.6000	1.6000
EDTA, Disodium Salt, Dihydrate	37.2500	37.2500
Ferrous Sulfate, Heptahydrate	27.8500	27.8500
Manganese Sulfate, Monohydrate	3.3300	3.3300
Potassium Iodide	0.8000	0.8000
Zinc Sulfate, Heptahydrate	1.5000	1.5000
Vitamins	mg/L	mg/L
Nicotinic Acid		0.5000
Pyridoxine, Hydrochloride		0.5000
Thiamine, Hydrochloride		1.0000