



Comparative Formulation

Earle's Balanced Salts

 Powdered Media	EBP01
Inorganic Salts	mg/L
Calcium Chloride, Dihydrate	264.9000
Magnesium Sulfate, Anhydrous	97.7000
Potassium Chloride	400.0000
Sodium Chloride	6800.0000
Sodium Phosphate, Monobasic, Monohydrate	140.0000
Other Components	mg/L
D-Glucose, Anhydrous	1000.0000
Phenol Red, Sodium Salt	17.0000



Liquid Media

EBL01

Inorganic Salts

mg/L

Calcium Chloride, Dihydrate	264.9000
Magnesium Sulfate, Anhydrous	97.7000
Potassium Chloride	400.0000
Sodium Bicarbonate	2200.0000
Sodium Chloride	6800.0000
Sodium Phosphate, Monobasic, Monohydrate	140.0000

Other Components

mg/L

D-Glucose, Anhydrous	1000.0000
Phenol Red, Sodium Salt	17.0000