

Dulbecco's PBS

Product Number PBL02-6X500ML

Dulbecco's phosphate-buffered saline, 1X. With calcium and magnesium. Store at room temperature.

Components	mg/L
Calcium Chloride, Anhydrous (CaCl ₂)	100.0000
Magnesium Chloride, Hexahydrate (MgCl ₂ · 6H ₂ O)	100.0000
Potassium Chloride (KCl)	200.0000
Potassium Phosphate, Monobasic, Anhydrous (KH ₂ PO ₄)	200.0000
Sodium Chloride (NaCl)	8000.0000
Sodium Phosphate Dibasic, Anhydrous (Na ₂ HPO ₄)	1150.0000