



CAISSON  
LABS

www.caissonlabs.com

## CERTIFICATE OF ANALYSIS

### Phosphate Buffered Saline

Catalog Number: PBP06  
Lot Number: 01240006  
MADE IN THE USA

Storage: 15° to 30°C  
Manufacture Date: Jan-2024  
Expiration Date: 31-Jan-2028

Test	Specifications	Results
Appearance	White Powder	Pass
Solubility in Water	9.57 g/L in ≤ 30 Minutes	7 Minutes
Moisture Content	≤ 2 %	0.4 %
Particle Size	≥ 90 % through 100 Mesh	83 %
pH ( <i>measured at 20°C</i> )	Record	QA: 7.13
Osmolality	Record mM/Kg	301 mM/Kg
Endotoxin	≤ 1.0 EU/mL	0.11 EU/mL
<i>Endotoxin Testing was Performed using the LAL Detection Method</i>		
Mycoplasma	Negative	Negative
<i>Mycoplasma Test was Performed using MycoAlert® Mycoplasma Detection Assay</i>		
Animal Cell Culture		
Toxicity	1 Pass	Pass
Functionality	≥ 80 %	93 %

**This product was produced and tested in accordance with all relevant SOP's and meets or exceeds the quality standards established by Caisson Laboratories, Inc.**

*This product is intended for laboratory use only.*

Prepared by:

Reviewed by:

Kurt Rieben  
QC Technician

26-Jan-2024

836 South 100 East  
Smithfield, Utah 84335

Phone: 435-755-7615  
FAX: 735-755-7617

BSE/TSE Declaration

To the present state of our knowledge we declare that this product does not contain any BSE/TSE AGENT. This declaration is given in good faith and no warranty express or implied with respect to quality and properties is made.

Toll Free: 1-877-840-0500