



CERTIFICATE OF ANALYSIS

Product Name: Phytoblend
Proprietary Blend of Agars

MADE IN THE USA

Catalog Number: PTP01

Storage: 15 to 30°C

Manufacture Date: MAR-2020

Lot Number: 03200014

Expiration Date: 31-MAR-2022

Test	Specifications	Results
Appearance	Off-White Powder	PASS
Visual Inspection	Product is Homogenous	PASS
Solubility in Water	9.00 g/L suspension	< 20 minutes
Moisture Content	≤ 22%	< 11 %
Particle Size (100 Mesh)	Record	100 %
pH (<i>measured at 26° C</i>)	Record	6.18
Gel Test	Gels at pH ≥ 5.5 after autoclaving and cooling	PASS
Plant Tissue Culture of: <i>Hosta spp.</i> <i>Solanum tuberosum</i> <i>Vitis spp.</i> <i>Tradescantia zebrina</i>	Supports plant growth with no morphological abnormalities	PASS

This product was produced and tested in accordance with all relevant SOP's and meets or exceeds the quality standards established by Caisson Laboratories, Inc.

This Product is intended for Laboratory Use Only.

Prepared by:

Zach Anthony
Quality Control

1-Apr-2020

Reviewed by:

BSE/TSE Declaration

To the present state of our knowledge we declare that this product does not contain any BSE/TSE AGENT. This declaration is given in good faith and no warranty express or implied with respect to quality and properties is made.