



www.caissonlabs.com

CERTIFICATE OF ANALYSIS

Phytoblend

Proprietary Blend of Agars. (Use at 7- 9 g/L)

Catalog Number: PTP01
Lot Number: 10200008
MADE IN THE USA

Storage: 15° to 30°C
Manufacture Date: Oct-2020
Expiration Date: **31-Oct-2022**

Test	Specifications	Results
Appearance	Off-White Powder	Pass
Solubility in Water	9.00 g/L in ≤ 30 minutes	6 minutes
Moisture Content	≤ 22 %	< 11 %
Particle Size	Record % through 100 Mesh	100 %
pH (<i>measured at 25°C</i>)	Record	7.21
Plant Tissue Culture of:		
<i>Hosta spp.</i>	Supports Plant Growth with	
<i>Solanum tuberosum</i>	No Morphological	Pass
<i>Vitis spp.</i>	Abnormalities	

This product was produced and tested in accordance with all relevant SOP's and meets or exceeds the quality standards established by Caisson Laboratories, Inc.

This product is intended for laboratory use only.

Prepared by:

Kurt Rieben
QC Technician

21-Oct-2020

Reviewed by:

BSE/TSE Declaration

To the present state of our knowledge we declare that this product does not contain any BSE/TSE AGENT. This declaration is given in good faith and no warranty express or implied with respect to quality and properties is made.