



Comparative Formulation

Woody Plant

 Powdered Media	WPP03	WPP04
Amino Acids	mg/L	mg/L
Glycine		2.0000
Macronutrients	mg/L	mg/L
Ammonium Nitrate	400.0000	400.0000
Calcium Chloride, Anhydrous	72.4700	72.4700
Calcium Nitrate, Tetrahydrate	556.0000	556.0000
Magnesium Sulfate, Anhydrous	180.7000	180.7000
Potassium Phosphate, Monobasic, Anhydrous	170.0000	170.0000
Potassium Sulfate, Anhydrous	990.0000	990.0000
Micronutrients	mg/L	mg/L
Boric Acid	6.2000	6.2000
Cupric Sulfate, Anhydrous	0.1600	0.1600
EDTA, Disodium Salt, Dihydrate	37.3000	37.3000
Ferrous Sulfate, Heptahydrate	27.8500	27.8500
Manganese Sulfate, Monohydrate	22.3000	22.3000
Molybdic Acid Sodium Salt, Dihydrate	0.2500	0.2500
Zinc Sulfate, Heptahydrate	8.6000	8.6000



Powdered Media

WPP03

WPP04

Vitamins

mg/L

mg/L

Myo-Inositol

100.0000

Nicotinic Acid

0.5000

Pyridoxine, Hydrochloride

0.5000

Thiamine, Hydrochloride

1.0000