CAISSON LABS
www.caissonlabs.com

GBP01 FORMULATION

Gamborg's B-5 Basal Salts

With Macronutrients and Micronutrients

Store at 2° to 8°C

Components	mg/L
Ammonium Sulfate $((NH_4)_2SO_4)$	134.00
Boric Acid (H_3BO_3)	3.00
Calcium Chloride, Anhydrous ($CaCl_2$)	113.24
Cobalt Chloride, Hexahydrate ($CoCl_2 \cdot 6H_2O$)	0.025
Cupric Sulfate, Pentahydrate ($CuSO_4 \cdot 5H_2O$)	0.025
EDTA, Disodium, Dihydrate ($C_{10}H_{14}N_2Na_2O_8 \cdot 2H_2O$)	37.26
Ferrous Sulfate, Heptahydrate ($FeSO4 \cdot 7H_2O$)	27.80
Magnesium Sulfate, Anhydrous ($MgSO_4$)	122.09
Manganese Sulfate, Monohydrate ($MnSO_4 \cdot H_2O$)	10.00
Molybdic Acid Sodium Salt, Dihydrate ($Na_2MoO_4 \cdot 2H_2O$)	0.250
Potassium Iodide (KI)	0.750
Potassium Nitrate (KNO_3)	2500.00
Sodium Phosphate, Monobasic, Monohydrate ($NaH_2PO_4 \cdot H_2O$)	150.08
Zinc Sulfate, Heptahydrate ($ZnSO_4 \cdot 7H_2O$)	2.00

836 South 100 East
Smithfield, Utah 84335

Phone: 435-755-7615
Email: custserv@caissonlabs.com

Toll Free: 1-877-840-0500