CAISSON LABS
www.caissonlabs.com

HBP01 FORMULATION

Hank's Balanced Salts

With Calcium, Magnesium and Phenol Red

Without Sodium Bicarbonate

Store at 15° to 30°C

Components	mg/L
Calcium Chloride, Anhydrous ($CaCl_2$)	140.00
D-Glucose, Anhydrous ($C_6H_{12}O_6$)	1000.00
Magnesium Sulfate, Anhydrous ($MgSO_4$)	98.00
Phenol Red, Sodium Salt ($C_{19}H_{13}NaO_5S$)	10.00
Potassium Chloride (KCl)	400.00
Potassium Phosphate Monobasic (KH_2PO_4)	60.00
Sodium Chloride (NaCl)	8000.00
Sodium Phosphate Dibasic, Anhydrous (Na_2HPO_4)	48.00

836 South 100 East
Smithfield, Utah 84335

Phone: 435-755-7615
Email: custserv@caissonlabs.com

Toll Free: 1-877-840-0500