www.caissonlabs.com

IVL01 FORMULATION

HEPES-TL

Supplement with BSA (Fract V), Sodium Pyruvate, and Gentamicin to make HEPES-TALP (Tyrode's Albumin Lactate Pyruvate). Store at 2° to 8°C.

Components	mg/L
Calcium Chloride, Anhydrous ($CaCl_2$)	220.20
HEPES ($C_8H_{18}N_2O_4S$)	2384.00
Magnesium Chloride, Hexahydrate ($MgCl_2 \cdot 6H_2O$)	83.90
Potassium Chloride (KCl)	238.72
Sodium Bicarbonate ($NaHCO_3$)	168.00
Sodium Chloride (NaCl)	6669.00
Sodium Phosphate, monobasic, Monohydrate ($NaH_2PO_4 \cdot H_2O$)	48.00
Sodium DL-Lactate, 60% ($C_3H_5NaO_3$) (Measured in mL)	1.416

836 South 100 East
Smithfield, Utah 84335

Phone: 435-755-7615
FAX: 735-755-7617

Toll Free: 1-877-840-0500