www.caissonlabs.com

PBL06 FORMULATION

1X Phosphate Buffered Saline Solution, pH 7.4

Without Calcium and Magnesium. Store at 15° to 30°C.

Components	mg/L
Potassium Phosphate Monobasic, Anhydrous (KH_2PO_4)	144.00
Sodium Chloride (NaCl)	9000.00
Sodium Phosphate Dibasic, Anhydrous (Na_2HPO_4)	421.62

836 South 100 East
Smithfield, Utah 84335

Phone: 435-755-7615
FAX: 735-755-7617

Toll Free: 1-877-840-0500